

# Investigation into factors affecting enzyme action

## Introduction

lodine is an indicator that turns blue/black when starch is present, but is otherwise brown. In this investigation a blue/black solution of starch and iodine will change to brown as the enzyme amylase digests/breaks down the starch into sugar.

The time taken for this reaction to occur is affected by temperature.

# **Apparatus**

test tube rack and six test tubes marker pen stopwatch 25 cm³ measuring cylinder 10 cm³ measuring cylinder beaker of 1% starch solution dropper bottle of iodine solution beaker of 10% amylase solution spotting tile dropping pipette

# Access to:

water bath or alternative method of heating water

## Method

- 1. Measure 10 cm<sup>3</sup> of 1 % starch solution into a test tube.
- 2. Measure 2 cm<sup>3</sup> of 10% amylase solution into a second test tube.
- 3. Place both tubes into a water bath set at 20 °C for 3 minutes.
- 4. Place a drop of iodine in six wells of a spotting tile.
- 5. Remove both test tubes from the water bath. Pour the amylase into the starch/iodine solution and start the stopwatch.
- 6. Immediately, use the dropping pipette to place one drop of the mixture onto the first drop of iodine. Record the colour of the solution.
- 7. Repeat step 6 every minute for five minutes.
- 8. Repeat steps 1-7 at 30°C, 40°C, 50°C, 60°C.

## **Analysis**

- 1. Use your observations to reach a conclusion regarding the effect of temperature on enzyme action.
- 2. Evaluate your method and suggest possible improvements.