

Simple extraction of DNA from living material

Introduction

DNA is the hereditary material found in all living things. In this practical you will extract the DNA from strawberries. Strawberries can have up to 8 copies of each chromosome and so contain a lot of DNA. When extracted from the strawberry this DNA is visible.

Apparatus

re-sealable plastic bag
 strawberry
 10cm³ washing up liquid (detergent)
 1 g sodium chloride
 100cm³ water
 2 × 250cm³ beakers (one beaker will be used for the filtering apparatus below)
 filter funnel
 coffee filter paper
 ice-cold 90% alcohol
 ice lolly stick or plastic coffee stirrer

Method

1. Remove the green top from the strawberry.
2. Put the strawberry into the plastic bag, seal it and crush for about 2 minutes.
3. In a beaker mix together 10cm³ of washing up liquid, 1 g of salt and 100cm³ water. This mixture is the DNA extraction liquid.
4. Add 10cm³ of the extraction liquid to the bag with the strawberry.
5. Re-seal the bag and gently mix the extraction liquid with the strawberry for 1 minute.
6. Place the coffee filter inside the beaker and gently pour the strawberry mixture into it.
7. Pour 10cm³ of ice-cold 90% ethanol down the side of the beaker into the strawberry mixture, do not mix or stir.
8. Within a few seconds you should see a white cloudy substance form in the clear layer above the strawberry mixture. Use a lolly stick to pull strands of this out of the top layer, this is the strawberry DNA.